




Product provides users of the nAb™ Agarose spin kit additional Spin Columns and Buffers. Quantities included are good for 20 reactions of nAb™ Agarose resin products. Use of spin columns eliminate inconsistent elutions that result from pipetting errors.


Fig 1 | Spin Column flow chart




Cell Lysate  
+  
nAb™  
Agarose



Binding to  
nAb™



Wash



agarose  
beads

Elute

## nAb™ Spin Kit Columns and Buffers,

spin columns and buffers for use with nAb™ Agarose Resin products.

### Box 1 | Product Info

nAb™ Spin Kit Columns and Buffers	
Cat.#: ABP-nAb-SPINK20	
Content	QTY
Spin Columns	20
Column top caps	20
Column bottom caps	40
Buffers	
10X Lysis Buffer	500µl
20X Binding Buffer	5ml
10X Wash Buffer	1.2ml
5X Elution Buffer	500µl
1X Neutralization Buffer	300µl

Store buffers at 4°C

### For Research Use Only. Not for Diagnostic or Therapeutic Use.

Purchase does not include or carry any right to re-sell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Allele Biotech is strictly prohibited

©2013 Allele Biotechnology and Pharmaceuticals, Inc. All rights reserved. The content mentioned herein are the property of Allele Biotechnology and Pharmaceuticals, Inc or their respective owners, except where otherwise stated. nAb™ is a registered trademark of Allele Biotechnology and Pharmaceuticals, Inc.