


# pmWasabi-CT-PuroR (ABP-FP-WCNCSP)

- CMV Promoter ..... 1-589
- Multiple Cloning Sites ..... 591-646
- mWasabi GFP ..... 653-1363
- SV40 poly(A) Signal ..... 1529-1579
- PGK Promoter ..... 1657-2159
- Puromycin Resistance Gene ..... 2242-2841
- SV40 poly(A) signal ..... 2965-3191
- Ampicillin Resistance Gene ..... 3841-4701
- pUC Origin ..... 4850-5492

Sequencing Primer:

CMV-F and mWasabi-R(ABP-FP-WASRPRM)



### Multiple Cloning Sites

```

541 ACGGTGGGAG GTCTATATAA GCAGAGCTGG TTTAGTGAAC CGTCAGATCC GCTAGCTCTA
601 GACTGGATCC TCGAATTCGG TACCTCAAGC TTGCGGCCGC CTCGAGGCCA CCATGGTGAG
661 CAAGGGCGAG GAGACCACAA TGGGCGTAAT CAAGCCCGAC ATGAAGATCA AGCTGAAGAT
    >>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>
1261 CATCGAGATC CTGAACCACG ACAAGGACTA CAACAAGGTG ACCGTTTACG AGATCGCCGT
1321 GGCCCGCAAC TCCACCGACG GCATGGACGA GCTGTACAAG TAA
  
```